




## Blisters

### Efficient long-term liquid reagent storage on lab-on-a-chip devices

Storage of liquid reagents on a microfluidic device is oftentimes a prerequisite for successful lab-on-a-chip applications. The use of blister pouches is not only a convenient way to incorporate liquids on-chip but also ensures long-term storage capability. Blister pouches serve as storage vessels and at the same time enable liquid movement through the application of pressure in a defined manner.



#### Off-The-Shelf Blister Pouches

- 50  $\mu$ l - 750  $\mu$ l filling volumes
- Compatible with inorganic solutions, acids and bases
- Reliable long-term liquid reagent storage
- Minimized contamination risk

#### Blister Test Chips

- Various off-the-shelf designs to evaluate blister emptying performance
- Feature *microfluidic ChipShop's* unique piercing blister seat technology for controlled liquid release

#### ChipGenie® edition BD


- High-precision blister driver instrumentation
- Compatible with most of *microfluidic ChipShop's* off-the-shelf blister test chips
- Controlled blister actuation operated with user-friendly control software and visual control unit

#### Customer-Specific Solutions & Services


- Blister filling service for customer-specific reagents
- Complex microfluidic cartridge development with comprehensive functional integration incl. liquid and dry reagent storage on-chip

## Blister Test Chip - Small Volume


microfluidic ChipShop offers a multitude of microfluidic chips to evaluate blister emptying. The blister test chip Fluidic 289 is particularly suited to investigate blisters with small volumes and comes with two blisters of 50  $\mu\text{l}$  and 100  $\mu\text{l}$  liquid content. For convenient visual assessment, both blisters are filled with dyed water.



Blister test chip with small volume blisters - blisters contain 50  $\mu\text{l}$  and 100  $\mu\text{l}$  of dyed water



Schematic drawing of small volume blister test chip Fluidic 289 with two blister seats and metering system - view from top



Schematic cross section of a standard blister setup on a microfluidic device


Product Code for Fluidic 289	Description	Blister details	Lid thickness [ $\mu\text{m}$ ]	Surface treatment	Material	Price [€/chip]		
						1+	10+	100+
10001305	Blister test chip Small volume	Blister position 1: 50 $\mu\text{l}$ ; cyan; $\text{H}_2\text{O}$ Blister position 2: 100 $\mu\text{l}$ ; magenta; $\text{H}_2\text{O}$	140	-	Topas	84.20	42.50	28.90
10001306	Blister test chip Small volume	Blister position 1: 50 $\mu\text{l}$ ; cyan; $\text{H}_2\text{O}$ Blister position 2: 100 $\mu\text{l}$ ; magenta; $\text{H}_2\text{O}$	140	hydrophilized	Topas	87.20	45.50	31.90

## Blister Test Chip - Volume Evaluation


Blister test chip Fluidic 522 comes with two mounted 250  $\mu\text{l}$  blisters, each of which is adjoined to a metering system. The blister emptying can be assessed with the aid of dyed blister content and either a chamber metering (blister position 1) or a channel metering system (blister position 2).




Blister test chip Fluidic 522 for volume evaluation - blisters contain 250  $\mu\text{l}$  of dyed water




Schematic drawing of blister test chip Fluidic 522 with two blister seats and metering systems - view from top



Schematic drawing of blister test chip Fluidic 522 - view from bottom (metering systems)



General components of a microfluidic chip exemplary shown by blister test chip Fluidic 522. The product consists of a chip body, attached blisters and a bonded cover lid



Close-up of blister seat with needles to pierce the bottom blister foil


Product Code for Fluidic 522	Description	Blister details	Lid thickness [ $\mu\text{m}$ ]	Surface treatment	Material	Price [€/chip]		
						1+	10+	100+
10001307	Blister test chip Volume evaluation	Blister position 1: 250 $\mu\text{l}$ ; cyan; $\text{H}_2\text{O}$ Blister position 2: 250 $\mu\text{l}$ ; magenta; $\text{H}_2\text{O}$	140	-	Topas	84.20	42.50	28.90
10001308	Blister test chip Volume evaluation	Blister position 1: 250 $\mu\text{l}$ ; cyan; $\text{H}_2\text{O}$ Blister position 2: 250 $\mu\text{l}$ ; magenta; $\text{H}_2\text{O}$	140	hydrophilized	Topas	87.20	45.50	31.90

Blister Test Chip - Volume Variation


This blister test chip is suitable for evaluating the emptying of a variety of blisters with volumes ranging from 50  $\mu$ l to 500  $\mu$ l. The chip design allows exact metering of 5  $\mu$ l (for 50  $\mu$ l blister), 15  $\mu$ l (for 150  $\mu$ l blister) and 25  $\mu$ l (for 250  $\mu$ l and 500  $\mu$ l blisters) with the aid of downstream channel/chambers.



Blister test chip Fluidic 1021 with blister pouches of various volumes. Each blister is associated to an appropriate metering system



Schematic drawing of blister test chip Fluidic 1021 with metering systems – view from top (blister seats and fluidic interfaces)



Schematic drawing of blister test chip Fluidic 1021 with metering systems – view from bottom (metering system)

Product Code for Fluidic 1021	Description	Blister details	Lid thickness [μm]	Surface treatment	Material	Price [€/chip]		
						1+	10+	100+
10001309	Blister test chip - Volume variation	Blister position 1: 500 μl; magenta; H <sub>2</sub> O Blister position 2: 250 μl; cyan; H <sub>2</sub> O Blister position 3: 150 μl; magenta; H <sub>2</sub> O Blister position 4: 50 μl; cyan; H <sub>2</sub> O	140	-	Topas	102.20	83.20	61.10
10001310	Blister test chip - Volume variation	Blister position 1: 500 μl; magenta; H <sub>2</sub> O Blister position 2: 250 μl; cyan; H <sub>2</sub> O Blister position 3: 150 μl; magenta; H <sub>2</sub> O Blister position 4: 50 μl; cyan; H <sub>2</sub> O	140	hydrophilized	Topas	105.20	86.20	64.10

Did you know?  
Controlled blister bursting


microfluidic ChipShop is using a unique blister seat technology. Thanks to piercing needles the blister content is released in a controlled and even manner. An erratic bursting of the blister is therefore obviated.

ChipGenie® edition BD


microfluidic ChipShop’s blister driver ChipGenie® edition BD is the tailored instrument for evaluation of blister emptying behavior. The device has been specifically developed to be used with microfluidic ChipShop’s microscope slide-format blister test chips, such as Fluidic 289 or Fluidic 522. The instrument ensures both a convenient, user-friendly chip insertion and a precise positioning of the inserted blister chip. Blister emptying is facilitated by two independent, high-precision vertical blister drives, while an integrated camera gives immediate visual feedback.

Instrument features:

- Two independent, high-precision vertical blister drives (driver speed 0.1  $\mu$ m/s up to 200  $\mu$ m/s; positioning accuracy +/- 20  $\mu$ m)
- Integrated camera for direct visual feedback
- Includes user-friendly software for easy system control
- Connection to PC established via USB port



Blisters of a blister test chip being actuated by the ChipGenie® edition BD instrument



Schematic succession of blister actuation by means of the ChipGenie® edition BD instrument

Product Code	Description	Dimensions [cm]	Price [€/instrument]
10000686	ChipGenie® edition BD – Blister Driver instrument	20 x 12 x 11	2.685,00




# Blisters - Liquid Reagent Storage


## Spare Blister Pouches

Blisters can be ordered as stand-alone parts being available with volumes ranging from 50  $\mu\text{l}$  to 750  $\mu\text{l}$  off-the-shelf. All blisters come with matching double-sided adhesive tape for convenient mounting on the microfluidic device. After use, blisters can be easily removed and replaced. All blisters can be ordered filled either with dye (cyan, magenta, yellow) or water (clear).

If you are interested in blisters filled with your customer-specific reagents please inquire our custom filling service. Simply write an e-mail with your requirements to [inquiries@microfluidic-ChipShop.com](mailto:inquiries@microfluidic-ChipShop.com).



Individual blister pouches with different filling volumes



Storage time dependent volume loss assessed for three different storage temperatures - blisters with 200  $\mu\text{l}$  ddH<sub>2</sub>O filling were stored at 4 °C, 20 °C and 37 °C and volumes were monitored over a time period of more than two years. At room temperature blister filling loss was less than 0.2% over 797 days

Product code (color-dependent)				Description	Blister Volume [ $\mu\text{l}$ ]	Price [€/blister]		
Cyan	Magenta	Yellow	Clear			1+	10+	100+
10001605	10001606	10001607	10001608	Blister pouch	50	8.20	5.80	2.95
10001610	10001611	10001612	10001613	Blister pouch	100	8.20	5.80	2.95
10001615	10001616	10001617	10001618	Blister pouch	150	8.20	5.80	2.95
10001620	10001621	10001622	10001623	Blister pouch	200	8.20	5.80	2.95
10001597	10001537	10001538	10001625	Blister pouch	250	8.20	5.80	2.95
10001602	10001601	10001600	10001599	Blister pouch	350	8.20	5.80	2.95
10001603	10001626	10001627	10001628	Blister pouch	500	8.20	5.80	2.95
10001665	10001666	10001667	10001668	Blister pouch	750	8.20	5.80	2.95


### Did you know? Liquid storage in tanks

Blisters are not the only way liquid reagents can be stored in microfluidics. In fact, for large volumes and short-term storage the use of tanks might be advisable. To learn more visit our webpage or have a look at our catalogue and our „Accessories“ leaflet!


## Tailored Blister Solutions

microfluidic ChipShop offers a variety of services on the subject of blister pouches. Looking to fill our standard blister pouches with reagents specific to your experiment? Simply inquire our blister filling service.

microfluidic ChipShop is also your specialized partner when it comes to the development of blister-containing microfluidic cartridges - from relatively simple customer-specific blister test chips to sophisticated, fully integrated lab-on-a-chip platforms. No matter how basic or complex your microfluidic assay requirements are, please contact us with your idea at [inquiries@microfluidic-ChipShop.com](mailto:inquiries@microfluidic-ChipShop.com)!



Simple blister test chip to evaluate blister emptying of a blister with a larger volume of 500  $\mu\text{l}$ . This test chip with Fluidic 761 is available as catalogue product with microfluidic ChipShop



Blister of different sizes mounted on a complex microfluidic cartridge, which has been developed at microfluidic ChipShop and features a great variety of functional units



A sophisticated lab-on-a-chip device with mounted blister, which are protected by a cover from unintentional actuation



FUJ1BLV3

